

## MOTION PRESERVATION MOTION PRESERVATION


PROCORAL TM









## **Anterior Cervical Disc Prosthesis**




SIZE	REF.CODE
12x14x5,0 mm	102.08 0250
12x14x5,5 mm	102.08 0255
12x14x6,0 mm	102.08 0260
12x14x6,5 mm	102.08 0265
12x14x7,0 mm	102.08 0270
12x14x7,5 mm	102.08 0275
12x14x8,0 mm	102.08 0280
12x14x8,5 mm	102.08 0285



REF.CODE
102.08 0150
102.08 0155
102.08 0160
102.08 0165
102.08 0170
102.08 0175
102.08 0180
102.08 0185



SIZE	REF.CODE
14x16x5,0 mm	102.08 0050
14x16x5,5 mm	102.08 0055
14x16x6,0 mm	102.08 0060
14x16x6,5 mm	102.08 0065
14x16x7,0 mm	102.08 0070
14x16x7,5 mm	102.08 0075
14x16x8,0 mm	102.08 0080
14x16x8,5 mm	102.08 0085



It is an alternative to the commonly performed anterior cervical discectomy and fusion (ACDF), a surgical procedure that is designed to address the pathology by eliminating motion at the diseased disc level.

Artificial disc replacement (ADR) surgery—also known as a total disc arthroplasty or total disc replacement (TDR)—is typically performed for a patient with a cervical disc herniation that is causing significant neck pain and/or arm pain that has not responded to nonsurgical treatment options and is significantly affecting the individual's quality of life and ability to function.

- Completely made of titanium and PEEK
- · PEEK material at the internal mechanism provides a super-smooth surface for a perfect motion capability
- PEEK cover which is shaped as a ring, placed around the internal mechanism in order to prevent the bone-fusion into the disc prosthesis
- One-piece and anatomical design
- Motion preservation
- Dark anodized and rough surface options are available
- Wide range of motion for flexion, extension, lateral bending and rotation
- Minimizing the scratching by the coating
- Smart design for the most effective movement capability
- Spikes for the enhanced hold on to vertebral bodies
- Reliable and stable connection with the instrument designed to withstand rotational forces and streamlined instrumentation
- Appropriate for Smith-Robinson approach
- Ti6Al4V (Grade 5) is used as raw materials for Prodorth Cervical Disc Prosthesis and these are originated only from our reliable partners in USA. PEEK material is supplied from EVONIK INDUSTRIES Germany. All certificates are available on request

